



Published in final edited form as:

Sci Transl Med. 2021 May 12; 13(593): . doi:10.1126/scitranslmed.abf6356.

A portable magnetofluidic platform for detecting sexually transmitted infections and antimicrobial susceptibility

Alexander Y. Trick¹, Johan H. Melendez², Fan-En Chen¹, Liben Chen³, Annet Onzia⁴, Aidah Zawedde⁴, Edith Nakku-Joloba⁵, Peter Kyambadde⁶, Emmanuel Mande⁴, Joshua Matovu⁴, Maxine Atuheirwe⁴, Richard Kwizera⁴, Elizabeth A. Gilliams^{2,7}, Yu-Hsiang Hsieh⁸, Charlotte A. Gaydos², Yukari C. Manabe^{2,4}, Matthew M. Hamill^{2,7}, Tza-Huei Wang^{1,2,3,9,*}

¹Department of Biomedical Engineering, Johns Hopkins University, Baltimore, MD 21218, USA.

²Division of Infectious Diseases, Johns Hopkins University School of Medicine, Baltimore, MD 21205, USA. ³Department of Mechanical Engineering, Johns Hopkins University, Baltimore, MD 21218, USA. ⁴Infectious Diseases Institute, Makerere University, Kampala, Uganda. ⁵School of Public Health, Makerere University, Kampala, Uganda. ⁶AIDS Control Program, Division of Sexually Transmitted Infections, Ministry of Health, Kampala, Uganda. ⁷Sexual Health Clinics, Baltimore City Health Department, Baltimore, MD 21205, USA. ⁸Department of Emergency Medicine, Johns Hopkins University School of Medicine, Baltimore, MD 21205, USA. ⁹Institute for NanoBiotechnology, Johns Hopkins University, Baltimore, MD 21218, USA.

Abstract

Effective treatment of sexually transmitted infections (STIs) is limited by diagnostics that cannot deliver results rapidly while the patient is still in the clinic. The gold standard methods for identification of STIs are nucleic acid amplification tests (NAATs), which are too expensive for widespread use and have lengthy turnaround times. To address the need for fast and affordable diagnostics, we have developed a portable, rapid, on-cartridge magnetofluidic purification and testing (PROMPT) polymerase chain reaction (PCR) test. We show that it can detect *Neisseria gonorrhoeae*, the pathogen causing gonorrhea, with simultaneous genotyping of the pathogen for resistance to the antimicrobial drug ciprofloxacin in <15 min. The duplex test was integrated into a

The Authors, some rights reserved; exclusive licensee American Association for the Advancement of Science. No claim to original U.S. Government Works

*Corresponding author. thwang@jhu.edu.

Author contributions:

T.-H.W. developed the project and guided the laboratory work. A.Y.T. and F.-E.C. designed the assay cartridge and instrument. A.Y.T. designed the mobile app and Bluetooth interface, and analyzed the data. J.H.M. and L.C. contributed to PCR assay design. A.Y.T. optimized the PCR assay for rapid thermocycling. A.Y.T. and J.H.M. conducted clinical sample testing. M.M.H., A.O., A.Z., and R.K. set up clinical study recruitment and survey administration. J.H.M. and E.M. cultured clinical samples and isolated bacterial strains for antimicrobial susceptibility testing. J.M. and M.A. processed clinical samples for nucleic acid testing. C.A.G. and E.A.G. designed the study in the BCHD clinics. Y.C.M., M.M.H., A.O., E.N.-J., and P.K. designed and coordinated the field study at the Infectious Diseases Institute in Kampala, Uganda. Y.-H.H., M.M.H., and T.-H.W. designed the end-user survey. A.Y.T. wrote the manuscript, and all co-authors edited the manuscript.

Competing interests:

A.Y.T. and T.-H.W. are co-inventors on patent PCT/US2019/029937 “A disposable reagent scaffold for biochemical process integration” that is associated with this study.

SUPPLEMENTARY MATERIALS

stm.sciencemag.org/cgi/content/full/13/593/eabf6356/DC1

[View/request a protocol for this paper from Bio-protocol.](#)

low-cost thermoplastic cartridge with automated processing of penile swab samples from patients using magnetic beads. A compact instrument conducted DNA extraction, PCR, and analysis of results while relaying data to the user via a smartphone app. This platform was tested on penile swab samples from sexual health clinics in Baltimore, MD, USA ($n = 66$) and Kampala, Uganda ($n = 151$) with an overall sensitivity and specificity of 97.7% (95% CI, 94.7 to 100%) and 97.6% (95% CI, 94.1 to 100%), respectively, for *N. gonorrhoeae* detection and 100% concordance with culture results for ciprofloxacin resistance. This study paves the way for delivering accessible PCR diagnostics for rapidly detecting STIs at the point of care, helping to guide treatment decisions and combat the rise of antimicrobial resistant pathogens.

INTRODUCTION

The rates of sexually transmitted infections (STIs) worldwide have increased markedly in the past decade, posing an enormous burden on global public health infrastructure (1). An estimated 1 million people worldwide acquire a curable STI daily (2). *Neisseria gonorrhoeae*, the etiological agent responsible for gonorrhea, is particularly concerning due to the rapid development of antimicrobial resistance (3, 4). With upward of 87 million annual *N. gonorrhoeae* infections worldwide, no available vaccine, a slow rate of antibiotic discovery, and a dearth of approved alternative antimicrobial treatments, the possibility of untreatable *N. gonorrhoeae* presages a looming global health care crisis (2, 4, 5).

Diagnosis of *N. gonorrhoeae* is the first step toward treatment and control. *N. gonorrhoeae* identification is routinely conducted using nucleic acid amplification tests (NAATs), such as polymerase chain reaction (PCR) (1, 6). Although NAATs have largely replaced culture-based methods as the gold standard for detection of *N. gonorrhoeae*, many NAATs take hours to complete, leading clinicians to prescribe antibiotic treatments empirically to avoid losing patients to follow-up. The cost of empiric treatment may be high, and a lack of diagnostic certainty can preclude use of more affordable oral medications for treating *N. gonorrhoeae* such as ciprofloxacin.

For low-income countries where testing options are absent or limited, particularly in sub-Saharan Africa, *N. gonorrhoeae* infections pose an especially heavy burden (2, 7). Limited STI surveillance and screening coupled with poor antimicrobial stewardship has led to a high prevalence of *N. gonorrhoeae* infections and antimicrobial resistance (7, 8). Rapid diagnostics at the point of care (POC) with antimicrobial resistance characterization have been highlighted as a viable solution to address this issue (9, 10). Fast, affordable, and accurate diagnostics that provide pathogen identification and antimicrobial resistance results are needed to shift the clinical workflow globally from syndromic treatment to precision therapy.

Commercialized rapid POC PCR technologies either are too slow, lack antimicrobial resistance data, or are too costly for widespread adoption in resource-limited settings (6, 11–16). We propose the use of magnetofluidic platforms coupled with inexpensive disposable assay cartridges to automate, accelerate, and simplify *N. gonorrhoeae* testing. Magnetofluidics leverages transport of functionalized magnetic beads to automate complex bioassays (17, 18). All manual sample handling steps typically required for a molecular

assay can be accomplished with magnetic field gradients using electromagnetic coils or by moving permanent magnets (19–21). We previously reported a magnetofluidic mobile NAAT platform with sample-to-answer detection of chlamydia in vaginal swabs, although the assay was limited in speed (>1 hour) and tested only a single genomic target for identification of the pathogen (22).

Here, we introduce the next generation of magnetofluidic cartridges and instrumentation and demonstrate the utility of our platform with a 15-min PCR assay for identifying *N. gonorrhoeae* and determining susceptibility of this pathogen to ciprofloxacin directly in clinical samples. This portable, rapid, on-cartridge magnetofluidic purification and testing (PROMPT) platform enabled automated on-site diagnosis in settings with confined workspaces, minimal resources, and unreliable electrical supply. We validated the PROMPT platform cartridge and compared it to culture assays and commercial NAAT assays using 35 clinical isolates and 66 penile swab eluates from sexual health clinics in Baltimore, MD, USA, as well as 151 penile swabs from multiple clinics in Kampala, Uganda. To assess the potential use and impact of the PROMPT platform, we performed on-site testing in sexual health clinics in Baltimore, conducted field testing in Kampala, and surveyed clinicians and medical laboratory technicians about POC testing in Kampala.

RESULTS

PROMPT platform design and operation

The PROMPT platform was designed for automated testing at the POC, which would allow clinicians to prescribe targeted antibiotic treatments, reduce overtreatment, and mitigate the need for follow-up visits (Fig. 1A). To enable POC use, the PROMPT platform was engineered to have a small footprint (12.7 cm by 13.4 cm by 8.4 cm) and require minimal power consumption. The instrument uses a portable mobile phone charger (5 V, ~2 A, 10,000 mAh), which supplies power to perform more than 20 tests.

The PROMPT workflow streamlined *N. gonorrhoeae* testing to three manual steps: elution of a swab, combining the swab eluate with a magnetic bead solution, and loading the combined solution into the cartridge (Fig. 1B). The disposable cartridge was the key feature of the platform that enabled automated sample preparation, rapid turnaround time, and a low cost per test. After sample loading, the cartridge was mounted onto the faceplate of the instrument (Fig. 1C and fig. S1). Within the cartridge, capture and release of the negatively charged nucleic acids was mediated by electrostatic attraction to a pH-responsive cationic coating on the magnetic beads (Fig. 1B). An acidic (pH 5) binding solution induced a positive charge on the beads for DNA capture and facilitated bacterial lysis (fig. S2). Transfer of the beads from the binding buffer into a neutral pH wash buffer removed contaminants that might inhibit the assay while retaining bound nucleic acids for subsequent elution directly into the PCR buffer (Fig. 1B). The more basic (pH 8.5) PCR buffer neutralized the beads' charge for release of the nucleic acids.

PCR-based detection of *N. gonorrhoeae* and determination of ciprofloxacin susceptibility were performed using an integrated two-color epifluorescence detector (23) for signal acquisition from duplexed hydrolysis probes targeting the multicopy *opa* gene (24) and

wild-type *gyrA* gene (25) (Fig. 1D). Fluorescence signals were wirelessly transferred to a smartphone app for analysis (Fig. 1D and movie S1). In the presence of an S91F mutation in the *gyrA* gene, which is highly predictive of ciprofloxacin resistance (26), the 6-carboxyfluorescein-labeled probe did not bind to the *gyrA* target sequence, resulting in absent signal amplification in the 6-carboxyfluorescein fluorescence channel. *N. gonorrhoeae*-positive samples with Cy5-labeled *opa* amplification and an absent *gyrA* signal were classified as ciprofloxacin resistant.

Magnetofluidic sample preparation and rapid PCR

The PROMPT platform achieved sample-to-answer diagnosis in <15 min through a combination of automated magnetofluidic sample preparation and rapid PCR. Transfer of the magnetic beads through the cartridge reagents was conducted with biaxial actuation of opposing neodymium magnets attached to a servo motor assembly (Fig. 2, A and B). Elution of the beads into the final reagent well coupled the sample preparation directly to the rapid PCR assay. The short turnaround time of the PCR assay was enabled by both enhanced assay chemistry to maximize reaction kinetics and optimized heat block geometry to quickly cycle the temperature.

Rapid PCR techniques overcome the necessity of long temperature hold times in traditional PCR routines by increasing the concentrations of primer, Mg^{2+} , and enzyme (DNA polymerase) to reduce the combined annealing and extension times to <1 s per cycle (27–29). Increased concentrations of primer, Mg^{2+} , and enzyme improved enzyme processivity and reduced the time needed for diffusion of reagents to anneal with target sequences such that primer extension occurred immediately upon reaching the primer annealing temperature (Fig. 2C). Denaturation of amplicons is a first-order process that occurs within 500 ms at sufficiently high temperatures (29). Previous studies on rapid PCR required specialized instruments impractical for clinical translation, overlooked sample preparation times, and used single-plex assays with intercalating dyes or electrophoresis for analysis (28–31). In contrast, our PROMPT cartridge integrated sample preparation with rapid PCR to form a multiplexed probe assay in a portable and user-friendly format. Compared to a traditional 40-cycle PCR assay with 20-s extension and 5-s denaturing hold times, the PROMPT assay's 1-s hold times translated to an overall decrease in annealing and denaturing time from 16.7 min to a little over 1 min (Fig. 2C).

By increasing the *opa* primer and enzyme concentrations to about seven- and fourfold of their respective standard concentrations (2 μ M primer, 1 U enzyme per 10- μ l reaction), we maximized the efficiency of the assay with 1-s hold times for annealing and denaturation (Fig. 2D). The specificity of the *opa* target amplification, verified by consistent melt curve peaks at 80°C (Fig. 2E), indicated specific amplification of the *opa* target, with increasing final amplicon quantity (taller melt curve peaks) correlating with increasing enzyme concentration. The 1 U per reaction was selected for the final assay as nonspecific amplification (melt peak aberration at 75°C) was noted with the 2.5 U concentration (Fig. 2E).

Even with ideal assay kinetics, a PCR assay is only as fast as the thermocycler supporting it. Traditional microfluidic devices require large and slow flat-bed heaters to match their planar

geometry (32), whereas our cartridge's extruded well design allowed for spatial isolation of the PCR reagents for targeted thermal control. To determine optimal geometries of the heating block for maximizing thermal ramp rates, finite element heat transfer simulations were conducted with a parametric sweep of heat block dimensions (figs. S3 and S4). These simulations demonstrated that at this scale, heat block volume was the primary determining factor in thermal ramp rate, and the thermal conductivity of aluminum ($\sim 200 \text{ W m}^{-1} \text{ K}^{-1}$) was sufficient to rapidly heat the PCR assay even with geometric constrictions resulting in cross-sectional areas as small as 0.05 mm^2 .

The final heat block design (Fig. 2F) minimized thermal mass by conforming precisely to the dimensions of the thermoelectric element with the well cantilevered to allow magnetic arm accessibility. A minimum heat block thickness of 0.25 mm was chosen due to practical machining limitations and to provide sufficient mechanical strength to withstand 5 lbs of force ($\sim 22.2 \text{ N}$) applied to the end of the well without yielding given a 10-fold safety margin. In the final assembly, heating rates ranged from 10.4° to 4.2°C/s and cooling rates from 10.0° to 18.5°C/s between 60° and 100°C (fig. S5), resulting in a <12 -min 40-cycle PCR compared to >50 -min 40-cycle PCR by conventional thermocyclers.

The use of elevated enzyme concentrations in a traditional PCR benchtop assay produced aberrations in the melt curves (black lines in Fig. 2G). However, with the PROMPT thermocycler, consistent melt curves with singular peaks were obtained regardless of enzyme concentration (blue lines in Fig. 2G). This demonstrated that the shorter annealing times used in the PROMPT PCR platform resulted in more specific amplification by limiting formation of less thermodynamically stable nonspecific dimers (27).

Assay cartridge design and validation

The cartridge wells contained the assay reagents in static droplets stabilized by surface tension using a continuous layer of silicone oil enabling captured analytes to be transferred via magnetic beads (Fig. 2B). Sample preparation within the cartridge was completed in 2.5 min, and 40 cycles of PCR were completed within 12 min, resulting in a <15 -min turnaround time (Fig. 3A). Traditional benchtop PCR thermocycler routines would complete <10 cycles in the same time period (Fig. 3B).

The cartridge design builds upon previous magnetofluidic cartridges (22, 23) assembled from laser-cut and thermoformed components with preloaded reagents (fig. S6). Combining the cost of assay reagents and thermoplastics resulted in a cost of US\$2.09 per cartridge (Fig. 3C). The use of polypropylene sheets for the thermoformed wells of the cartridge provided a chemically inert PCR-compatible scaffold, which allowed for 1 month of storage with consistent performance (fig. S7 and data file S1).

The limit of detection (LoD) of the cartridge assay was evaluated using serial dilutions of a ciprofloxacin-susceptible *N. gonorrhoeae* reference strain [American Type Culture Collection (ATCC) 49226]. The *opa* assay for detection of *N. gonorrhoeae* showed successful amplification of *opa* in all replicates down to 10 colony-forming units (CFUs); the *gyrA* assay for detection of ciprofloxacin resistance demonstrated a less sensitive LoD of 100 CFUs (Fig. 3, D and E). Real-time reporting of results using a custom regression

algorithm (figs. S8 and S9) allowed for a turnaround time of <10 min for samples with a high bacterial load (Fig. 3F).

The performance of the *opa* and *gyrA* assays was evaluated against a panel of well-characterized *N. gonorrhoeae* reference strains and archived *N. gonorrhoeae* clinical isolates (Fig. 3G and table S1). Only one of the *N. gonorrhoeae* isolates produced a false-negative *opa* result (BCHD-23); the remaining 37/38 *N. gonorrhoeae* isolates were identified by the *opa* assay. All ciprofloxacin-susceptible isolates produced a detectable wild-type *gyrA* amplification signal.

Specificity of the cartridge assay was assessed with 15 bacterial and viral species known to infect the human urogenital tract, including six *Neisseria* species and four *Chlamydia trachomatis* serovars (Fig. 3G and table S2). None of the pathogens produced any false-positive amplification with the *opa* assay. Cross-reactivity of the *gyrA* assay was observed with the bacterial species *Neisseria mucosa*. The *gyrA* gene may also be present in other bacteria, particularly in *Neisseria* commensal species. In order for positive *gyrA* amplification to be considered valid for detection of ciprofloxacin-susceptible *N. gonorrhoeae* strains, the PROMPT assay must also be positive for *opa* gene amplification.

Validation of the PROMPT assay using clinical samples

Before blinded testing or field testing, the clinical sensitivity of the PROMPT cartridge assay was evaluated using 32 archived urethral swab eluates (33). All samples were correctly classified as *N. gonorrhoeae* positive ($n = 18$) or *N. gonorrhoeae* negative ($n = 14$) by the PROMPT assay in full concordance with comparator results from NAAT assays (table S3). Of the 18 *N. gonorrhoeae*-positive samples, 13 were identified as having *gyrA* wild-type sequences indicative of ciprofloxacin susceptibility. Phenotypic susceptibility data were not available for these samples.

After platform validation with archived clinical samples, field testing at the POC using a blinded approach was implemented at two Baltimore City Health Department (BCHD) sexual health clinics during May and June 2019. Penile-meatal clinical swab samples were prospectively collected after informed consent, and tested on the PROMPT platform. Of the 34 samples tested, two were identified as *N. gonorrhoeae* positive by the PROMPT cartridge assay, by culture, and by commercial NAAT assays. All *N. gonorrhoeae*-negative samples were correctly classified by the PROMPT cartridge assay. The two *N. gonorrhoeae*-positive samples were identified as having wild-type *gyrA* sequences. Phenotypic antimicrobial susceptibility testing confirmed that the isolates cultured from these two samples were susceptible to ciprofloxacin. Figure 4A illustrates the PROMPT cartridge sensitivity with the *opa* assay for *N. gonorrhoeae* detection using the Baltimore clinic swabs (orange and light gray circles). Figure 4B shows the error-free prediction of ciprofloxacin susceptibility with the *gyrA* assay in the Baltimore clinic strains isolated from swabs and based on antimicrobial resistance profiles determined by culture.

Additional blinded testing of the PROMPT cartridge assay was performed using participant-collected penile-meatal swabs after informed consent, from six clinics in Kampala, Uganda between October 2019 and February 2020 ($n = 151$) (fig. S10). Swab eluates ($n =$

51) were tested on-site at the Infectious Diseases Institute (IDI) in Kampala with two PROMPT instruments and compared to culture and commercial NAAT results. Because of the coronavirus disease 2019 (COVID-19) pandemic, the remaining swabs ($n = 100$) were analyzed in the laboratory at Johns Hopkins University. Of the 151 samples collected in Kampala, 114 were *N. gonorrhoeae* positive by a NAAT assay. Two samples were negative by NAAT but confirmed as *N. gonorrhoeae* positive by culture, resulting in a prevalence of 76.8% (116 of 151). The PROMPT cartridge assay failed to detect three *N. gonorrhoeae*-positive samples, and two *N. gonorrhoeae*-negative samples were classified as false positives (Fig. 4A). All cultured *N. gonorrhoeae* isolates from the Kampala samples ($n = 85$) were determined to be ciprofloxacin resistant by antimicrobial susceptibility testing, which is consistent with published local antimicrobial resistance patterns (7), and in full concordance with results of the *gyrA* assay.

Receiver operating characteristic (ROC) curves for the identification of *N. gonorrhoeae* by *opa* genotyping and antimicrobial susceptibility testing by *gyrA* genotyping (Fig. 4, C and D) had area under the curve values of 0.987 [95% confidence interval (CI), 0.973 to 1.000] and 1.0, respectively. The PROMPT cartridge assay matched the NAAT and culture results with a high degree of concordance for an overall clinical sensitivity and specificity of 97.7% (95% CI, 94.7 to 100%) and 97.6% (95% CI, 94.1 to 100%), respectively. Results obtained using the PROMPT cartridge assay were compared to NAAT, culture, and phenotypic antimicrobial susceptibility testing results. (See data file S2 for compiled results from all cartridge runs and table S4 for a summary of discrepant results.) The average full assay turnaround time was 14.5 min (SD, 0.9 min; median, 14.4 min), starting with the cartridge mounted into the PROMPT instrument to completion of 40 PCR cycles (Fig. 4E). When using the live-reporting algorithm to determine when amplification was first evident, the turnaround time decreased to 11.1 min on average (SD, 1.5 min; median, 10.9 min).

End-user survey

To assess the potential utility of POC testing with the PROMPT platform, we conducted a survey of clinicians and laboratory staff ($n = 30$) in Kampala (figs. S11 and S12 and data file S3). Most (93%) of the respondents indicated that assay sensitivity was the most important characteristic of a diagnostic test for *N. gonorrhoeae*, followed by turnaround time (76.6%). For POC testing in general, participants preferred a test that was fast (24 of 30) and simple to use (22 of 30), with less than 50% (13 of 30) prioritizing cost or mobility of the test (12 of 30) (fig. S11B).

The survey also assessed general attitudes about POC testing, duration of typical patient visits, and desired test duration to fit clinical workflow (fig. S11C). POC tests were generally preferred to laboratory testing (24 of 30), and more than 50% (16 of 30) of respondents stated that a typical patient's visit lasted <1 hour. To fit POC testing into the timeline of a patient's visit, most respondents desired a turnaround time of <30 min (23 of 30); six participants preferred a turnaround time of <15 min.

After indicating preferences for desired testing characteristics, participants were shown a demonstration of the PROMPT platform operation (movie S1), after which they answered questions regarding their impressions of the test. At least 90% of survey participants

selected “extremely easy” or “easy” when asked about ease of following instructions, running the test, and interpreting results. Participants commented positively that such a test would save them from “long waits for the cumbersome gram [sic] stain and culture tests”; resistance profiles were highly desired to replace their current tests “which may take...almost a week in our health centre to get results.” One respondent was concerned about the utility of ciprofloxacin resistance testing in Uganda due to the high prevalence of antimicrobial resistance, and suggested cephalosporin susceptibility testing instead. Other concerns pertained to the limited availability of smartphones, and that connection to a computer would improve patient data security and the ability to retrieve results. Despite the few concerns raised, nearly all of the participants (29 of 30) indicated that they “Definitely” or “Probably” would use the PROMPT platform if it was made available.

DISCUSSION

Rapid identification of the bacterial pathogen *N. gonorrhoeae* and associated antimicrobial susceptibility will help to promote better antimicrobial stewardship to reduce selection of antimicrobial resistant bacterial strains. Although ciprofloxacin is no longer recommended for the treatment of *N. gonorrhoeae* in the United States (34), up to 80% of *N. gonorrhoeae* infections are susceptible to ciprofloxacin (35). POC use of our PROMPT platform to rapidly identify ciprofloxacin-susceptible *N. gonorrhoeae* infections opens up opportunities for providing access to less expensive oral ciprofloxacin compared to the current first-line recommended injectable antibiotic ceftriaxone (36–39).

Recent POC diagnostics have trended toward isothermal techniques, whereas our work demonstrates that PCR-based assays can provide a fast turnaround time while maintaining high sensitivity and specificity in a cost-effective POC format. A major advantage of using a proven NAAT strategy like PCR over isothermal techniques is the opportunity for facile adoption of established clinical assays with the PROMPT cartridge. The *opa* and *gyrA* assays that we used in our PROMPT cartridge have been shown to be highly sensitive and specific, with <1% false-positive or false-negative results (24, 40–45).

Proponents of isothermal NAAT platforms as alternatives to PCR highlight the presumed advantages in instrument simplicity or “instrument-free” operation, although many fail to emphasize the necessity for nucleic acid purification before testing (46–48). Without integrated sample preparation, these techniques require a trained professional and laboratory equipment (heaters, centrifuges, etc.) to process clinical samples, which ultimately provides little benefit over traditional PCR assays. Our PROMPT platform addresses both sample preparation and rapid turnaround time in a format that requires minimal training to operate. Tests with instrumentation and connectivity capabilities may be preferable in clinical settings to remove risk of errors associated with user interpretation of test readouts or for interfacing with health care records (as noted by our survey respondents). By avoiding complex moving parts and precision labyrinthine microchannels traditionally found in POC diagnostic devices, the magnetofluidic cartridge is suitable for high-throughput manufacture at an affordable cost. Compared to the leading industry POC PCR diagnostic cartridges, our magnetofluidic cartridge design has potential to provide even faster results, contribute less

plastic waste to the environment, and enable greater access to diagnostics in low-resource settings.

For translation of the PROMPT platform to patient care settings, there are limitations in our study that must be addressed. Performance of the test with vaginal and extragenital swabs, and urine samples needs to be validated. Future cartridges would benefit from higher multiplexing to include assay controls, other antimicrobial resistance markers, or additional STIs such as chlamydia. User workflow could be further simplified by incorporating swab elution into the magnetic bead mixture directly in the cartridge. Furthermore, cartridges in the current format are limited to cold chain storage (-20°C). For widespread adoption, it may be necessary to adapt the cartridges for dry reagents to ensure prolonged shelf-life at room temperature. Nonetheless, advances demonstrated so far in magnetofluidic cartridge technology show promise for delivering rapid gold standard infectious disease diagnostics in a scalable and accessible format that could improve clinical management of STIs.

MATERIALS AND METHODS

Study design

We performed a cross-sectional, nonrandomized diagnostic accuracy study to determine the analytical performance of our PROMPT platform for detection of *N. gonorrhoeae* and characterization of ciprofloxacin susceptibility, and compared its performance with that of laboratory-based NAAT and culture assays. Following test development and preliminary testing using archived specimens, male participants were offered enrollment to self-collect penile-meatal swabs to be tested on the PROMPT platform and in comparator gold standard tests. Following written informed consent, male participants were recruited in Baltimore, MD, USA and Kampala, Uganda. Whenever possible, testing of the new assay was performed on-site to evaluate the performance in clinical settings outside of the research laboratory. PROMPT operators were blinded to NAAT and culture results. A total of 34 participants in Baltimore and 250 participants in Kampala were recruited and provided samples for the study, but only 151 samples from Kampala were included in this study due to delays associated with the ongoing COVID-19 pandemic. Ethical oversight was provided by the Johns Hopkins Institutional Review Board (IRB) for both studies (IRB numbers 00215298 and 00199514). In addition, the Joint Clinical Research Center at Makerere University (protocol reference number JC0919) and the Ugandan National Council for Science and Technology (study number HS455ES) approved the study.

For the Baltimore study, we enrolled a prospective, convenience sample of participants to evaluate the performance of the PROMPT platform in a clinical setting with low prevalence of *N. gonorrhoeae* and ciprofloxacin resistance. This initial study was used to develop the protocol for testing prospectively collected samples in a clinical setting before the Uganda study. For the Uganda study, to obtain the minimal sample size, the area under the ROC curve (AUC) with the PROMPT assay was estimated. In these estimations, *N. gonorrhoeae* positivity in men presenting to clinics in Kampala with urethral discharge was assumed to be between 75 and 85%. Using a null hypothesis of AUC at 0.80, assuming a 0.10 increase in AUC with the PROMPT assay and with a type 1 error of 5% and 80% power, a total of 213 participants was required. However, anticipating that 15% or greater results would be

uninterpretable (the test result was invalid or indeterminate), the sample size was adjusted to a total of 250 study participants.

Assay reagents

Magnetic particle binding buffer solutions for each cartridge were formulated using 10 μl of acetic acid-based binding buffer, 4 μl of ChargeSwitch magnetic particles (25 mg/ml) (ChargeSwitch gDNA Mini Bacteria Kit, Invitrogen), and 1 or 10 μl of 10% Tween 20 (fig. S13). For magnetic bead wash buffer, the wash buffer from the same ChargeSwitch kit was supplemented with a final concentration of 0.2% (v/v) Tween 20 to reduce surface tension between the buffer and the oil and permit maximum magnetic bead extraction. The PCR buffers and sequences of the primers and probes are described in tables S5 and S6. Direct elution of the magnetic beads into the PCR buffers was shown to remove primers and enzyme from the assay (fig. S14), but these effects were compensated with excess primer and enzyme in the PROMPT PCR.

PROMPT cartridge

The magnetofluidic cartridges were composed of three thermoplastic layers (fig. S6). The bottom layer was fabricated by thermoforming 0.2-mm-thick polypropylene sheet (AKAHA) over three-dimensional (3D)-printed molds (Form 2, USA) designed in Solidworks 2017 computer-aided design (CAD) software (Dassault Systèmes) to produce extruded wells. The middle layer was laser-cut from 1.5-mm-thick clear acrylic [polymethyl methacrylate (PMMA)] sheet (McMaster-Carr) with pressure-sensitive adhesive (PSA) (9472LE adhesive transfer tape, 3M) laminated on both sides. The top layer was laser-cut from 0.75-mm-thick acrylic (ePlastics) with polytetrafluoroethylene (PTFE) tape (McMaster-Carr) laminated to one side and patterned by laser-etching.

Assay reagents (10 μl of PCR solution and 50 μl of wash buffer) were loaded into the corresponding wells of the thermoformed section, followed by joining the thermoformed section with PSA, and loading the cartridge with 420 μl of silicone oil (50 c St, Millipore-Sigma) to cover the wells and remaining empty space. The cartridges were used immediately, or the sample injection port was sealed with adhesive tape (Scotch Magic Tape, 3M) and the cartridge was stored at -20°C until use. For testing, an aliquot of sample and magnetic bead binding buffer (acceptable volumes between 10 and 100 μl) was injected directly into the first well of the cartridge and the cartridge was mounted onto the PROMPT instrument faceplate with the PCR well (Fig. 1C).

Instrumentation design and assembly

The housing of the PROMPT instrument and magnetic arm was 3D-printed (Formlabs Form 2, black resin) and equipped with permanent neodymium magnets for magnetic clamping, alignment, and magnetofluidic transfer in the cartridges (fig. S1). A microcontroller (Arduino Uno R3) monitored and controlled the heat block temperature with current supplied by a motor shield (Arduino Motor Shield Rev3), and an additional custom printed circuit board shield provided wiring connections for the fluorescence detector (Fluo Sens Integrated, Qiagen) and magnetic arm motors. The magnetic arm motors consisted of a linear servo (PQ12-R, Actuonix), which was secured to the shaft of a rotary servo

(HS-485HB Hitec RCD) for biaxial actuation. The fluorescence detection and heating module was assembled as previously described (23), with the exception of the design of the aluminum heat block that was optimized for rapid thermocycling. Finite element heat transfer simulations to characterize effects of heat block dimensions on temperature ramp rates were conducted using COMSOL Multiphysics software version 5.1.

The instrument communicated to our custom smartphone app through a wireless serial connection using a Bluetooth module (Bestgle HC-05) connected to the microcontroller. A switch in the back of the instrument controlled power to the Bluetooth module to switch the communication from smartphone/tablet mode to universal serial bus (USB) connection for programming with a computer. The instrument was powered through a barrel-plug connection to a portable 5-V power bank (Anker PowerCore, 10,000 mAh) or an alternating current (AC) adapter.

To control the instrument and automate data collection and analysis, a graphic user interface (GUI) was developed in Java using Processing (Processing 3.4, www.processing.org) for both PC and Android devices. During clinical testing, the devices were controlled with a Samsung Galaxy S4 smartphone or Android tablet (Asus ZenPad). The GUI displayed real-time fluorescence data and live updates if amplification for either target was detected (movie S1).

PCR real-time fluorescence analysis algorithm

A custom algorithm for identifying the presence of PCR amplification was implemented in the Processing GUI to automatically report the *opa* and *gyrA* results in real time. The algorithm is described in figs. S8 and S9.

Clinical isolates and archived swab testing

N. gonorrhoeae isolates from international repositories (ATCC 49226, WHO K, WHO L) or obtained from previous studies in the BCHD clinics ($n = 35$) (49) were cultured as previously described (48) and tested on the PROMPT cartridge to evaluate assay sensitivity. Phenotypic ciprofloxacin susceptibility data were available for all *N. gonorrhoeae* isolates (table S1). To evaluate cartridge functionality in a clinical matrix, archived eluates of well-characterized urethral swabs collected from BCHD clinics in 2015 (33) were also tested on the PROMPT platform.

Clinical enrollment and sample collection

After providing written informed consent, patients in the Baltimore, MD, study self-collected two penile-swabs; one of the swabs was eluted into 500 μ l of phosphate-buffered saline (PBS). The PBS eluate was used for PROMPT analysis (50 μ l), *N. gonorrhoeae* culture (100 μ l), and NAAT assay (Hologic Aptima CT/NG) (200 μ l), and the remainder was stored at -20°C for retesting or subsequent studies.

In the Kampala, Uganda study, after informed consent, participants self-collected two penile-meatal swabs at six local health clinics. Before study enrollment, participants provided a urethral swab for the World Health Organization (WHO) Enhanced Gonococcal

Antimicrobial Surveillance Programme (EGASP) (fig. S10), which was cultured for *N. gonorrhoeae* and tested for antimicrobial susceptibility. At the Makerere University IDI laboratory, the penile-meatal swabs were each eluted into 500 µl of PBS, and the eluates were combined and used for PROMPT testing and NAAT. Samples were either tested immediately on PROMPT cartridges or stored frozen until tested.

Statistical analysis

Data analysis of clinical results was conducted using R software version 3.6.2. Generation of ROC curves with sensitivity, specificity, and AUC calculations was done using the pROC package version 1.16.2.

Supplementary Material

Refer to Web version on PubMed Central for supplementary material.

Acknowledgments:

We thank M. Barnes and P. Barnes for help in setting up our study and on-site testing in the BCHD clinics. We also thank T. Hu and J. Chung for assistance in cartridge and instrument testing.

Funding:

This work was funded by the NIH (grants R01AI138978, U54EB007958, and R61AI154628 to T.-H.W.).

Data and materials availability:

The software for the mobile application programmed in Processing (version 3.4) including the algorithm for fluorescence signal analysis is available at <https://doi.org/10.5281/zenodo.4613082>. All other data associated with this study are present in the paper or the Supplementary Materials.

REFERENCES AND NOTES

1. Hull S, Kelley S, Clarke JL, Sexually transmitted infections: Compelling case for an improved screening strategy. *Popul. Health Manag*20, S-1–S-11 (2017).
2. Rowley J, Vander Hoorn S, Korenromp E, Low N, Unemo M, Abu-Raddad LJ, Chico RM, Smolak A, Newman L, Gottlieb S, Thwin SS, Broutet N, Taylor MM, Chlamydia, gonorrhoea, trichomoniasis and syphilis: Global prevalence and incidence estimates, 2016. *Bull. World Health Organ*97, 548–562 (2019). [PubMed: 31384073]
3. Unemo M, Golparian D, Eyre DW, Antimicrobial resistance in *Neisseria gonorrhoeae* and treatment of gonorrhea, in *Neisseria gonorrhoeae*, *Methods in Molecular Biology*, Christodoulides M, Ed. (Humana, 2019), vol. 1997, pp. 37–58. [PubMed: 31119616]
4. Unemo M, Shafer WM, Antimicrobial resistance in *Neisseria gonorrhoeae* in the 21st Century: Past, evolution, and future. *Clin. Microbiol. Rev*27, 587–613 (2014). [PubMed: 24982323]
5. Alirol E, Wi TE, Bala M, Bazzo ML, Chen XS, Deal C, Dillon JAR, Kularatne R, Heim J, Hooft van Huijsduijnen R, Hook EW, Lahra MM, Lewis DA, Ndowa F, Shafer WM, Tayler L, Workowski K, Unemo M, Balasegaram M, Multidrug-resistant gonorrhea: A research and development roadmap to discover new medicines. *PLOS Med.* 14, e1002366 (2017). [PubMed: 28746372]
6. Gaydos CA, Melendez JH, Point-by-point progress: Gonorrhea point of care tests. *Expert Rev. Mol. Diagn*20, 803–813 (2020). [PubMed: 32500764]

7. Workneh M, Hamill MM, Kakooza F, Mande E, Wagner J, Mbabazi O, Mugasha R, Kajumbula H, Walwema R, Zenilman J, Musinguzi P, Kyambadde P, Lamorde M, Manabe YC, Antimicrobial resistance of *Neisseria gonorrhoeae* in a newly implemented surveillance program in Uganda: Surveillance report. *JMIR Public Health Surveill.* 6, e17009 (2020). [PubMed: 32519969]
8. Fingerhuth SM, Bonhoeffer S, Low N, Althaus CL, Antibiotic-resistant *Neisseria gonorrhoeae* spread faster with more treatment, not more sexual partners. *PLOS Pathog.* 12, e1005611 (2016). [PubMed: 27196299]
9. Tuite AR, Gift TL, Chesson HW, Hsu K, Salomon JA, Grad YH, Impact of rapid susceptibility testing and antibiotic selection strategy on the emergence and spread of antibiotic resistance in gonorrhea. *J. Infect. Dis*216, 1141–1149 (2017). [PubMed: 28968710]
10. Sadiq ST, Mazzaferri F, Unemo M, Rapid accurate point-of-care tests combining diagnostics and antimicrobial resistance prediction for *Neisseria gonorrhoeae* and *Mycoplasma genitalium*. *Sex. Transm. Infect*93, S65–S68 (2017). [PubMed: 28684610]
11. Gaydos CA, Review of use of a new rapid real-time PCR, the Cepheid GeneXpert® (Xpert) CT/NG assay, for *Chlamydia trachomatis* and *Neisseria gonorrhoeae*: Results for patients while in a clinical setting. *Expert Rev. Mol. Diagn*14, 135–137 (2014). [PubMed: 24450867]
12. Causer LM, Guy RJ, Tabrizi SN, Whiley DM, Speers DJ, Ward J, Tangey A, Badman SG, Hengel B, Natoli LJ, Anderson DA, Wand H, Wilson D, Regan DG, Shephard M, Donovan B, Fairley CK, Kaldor JM, Molecular test for chlamydia and gonorrhoea used at point of care in remote primary healthcare settings: A diagnostic test evaluation. *Sex. Transm. Infect*94, 340–345 (2018). [PubMed: 29748180]
13. Widdice LE, Hsieh Y-HH, Silver B, Barnes M, Barnes P, Gaydos CA, Performance of the atlas genetics rapid test for *Chlamydia trachomatis* and women's attitudes toward point-of-care testing. *Sex. Transm. Dis*45, 723–727 (2018). [PubMed: 29771869]
14. Low N, Unemo M, Molecular tests for the detection of antimicrobial resistant *Neisseria gonorrhoeae*: When, where, and how to use? *Curr. Opin. Infect. Dis*29, 45–51 (2016). [PubMed: 26658656]
15. Unemo M, Bradshaw CS, Hocking JS, de Vries HJC, Francis SC, Mabey D, Marrazzo JM, Sonder GJB, Schwebke JR, Hoornenborg E, Peeling RW, Philip SS, Low N, Fairley CK, Sexually transmitted infections: Challenges ahead. *Lancet Infect. Dis*17, e235–e279 (2017). [PubMed: 28701272]
16. Gaydos C, Hardick J, Point of care diagnostics for sexually transmitted infections: Perspectives and advances. *Expert Rev. Anti Infect. Ther*12, 657–672 (2014). [PubMed: 24484215]
17. Shin DJ, Wang TH, Magnetic droplet manipulation platforms for nucleic acid detection at the point of care. *Ann. Biomed. Eng*42, 2289–2302 (2014). [PubMed: 25008142]
18. Zhang Y, Nguyen N-T, Magnetic digital microfluidics—A review. *Lab Chip*17, 994–1008 (2017). [PubMed: 28220916]
19. Lehmann U, Vandevyver C, Parashar VK, Gijs MAM, Droplet-based DNA purification in a magnetic lab-on-a-chip. *Angew. Chem. Int. Ed. Engl*45, 3062–3067 (2006). [PubMed: 16555359]
20. Zhang Y, Wang TH, Full-range magnetic manipulation of droplets via surface energy traps enables complex bioassays. *Adv. Mater*25, 2903–2908 (2013). [PubMed: 23529938]
21. Chiou C-H, Jin Shin D, Zhang Y, Wang T-H, Topography-assisted electromagnetic platform for blood-to-PCR in a droplet. *Biosens. Bioelectron*50, 91–99 (2013). [PubMed: 23835223]
22. Shin DJ, Athamanolap P, Chen L, Hardick J, Lewis M, Hsieh YH, Rothman RE, Gaydos CA, Wang TH, Mobile nucleic acid amplification testing (mobiNAAT) for *Chlamydia trachomatis* screening in hospital emergency department settings. *Sci. Rep*7, 4495 (2017). [PubMed: 28674410]
23. Shin DJ, Trick AY, Hsieh Y-H, Thomas DL, Wang T-H, Sample-to-answer droplet magnetofluidic platform for point-of-care hepatitis C viral load quantitation. *Sci. Rep*8, 9793 (2018). [PubMed: 29955160]
24. Tabrizi SN, Chen S, Tapsall J, Garland SM, Evaluation of opa-based real-time PCR for detection of *Neisseria gonorrhoeae*. *Sex. Transm. Dis*32, 199–202 (2005). [PubMed: 15729160]

25. Giles J, Hardick J, Yuenger J, Dan M, Reich K, Zenilman J, Use of applied biosystems 7900HT sequence detection system and taqman assay for detection of quinolone-resistant *Neisseria gonorrhoeae*. *J. Clin. Microbiol*42, 3281–3283 (2004). [PubMed: 15243093]
26. Allan-Blitz LT, Wang X, Klausner JD, Wild-type gyrase a genotype of *Neisseria gonorrhoeae* predicts in vitro susceptibility to ciprofloxacin: A systematic review of the literature and meta-analysis. *Sex. Transm. Dis*44, 261–265 (2017). [PubMed: 28407640]
27. Wittwer CT, Garling DJ, Rapid cycle DNA amplification: Time and temperature optimization. *Biotechniques*10, 76–83 (1991). [PubMed: 2003928]
28. Farrar JS, Wittwer CT, Extreme PCR: Efficient and specific DNA amplification in 15–60 seconds. *Clin. Chem*61, 145–153 (2015). [PubMed: 25320377]
29. Millington AL, Houskeeper JA, Quackenbush JF, Trauba JM, Wittwer CT, The kinetic requirements of extreme qPCR. *Biomol. Detect. Quantif*17, 100081 (2019). [PubMed: 31285997]
30. Neuzil P, Zhang C, Pipper J, Oh S, Zhuo L, Ultra fast miniaturized real-time PCR: 40 cycles in less than six minutes. *Nucleic Acids Res.* 34, e77 (2006). [PubMed: 16807313]
31. Wheeler EK, Hara CA, Frank J, Deotte J, Hall SB, Benett W, Spadaccini C, Beer NR, Under-three minute PCR: Probing the limits of fast amplification. *Analyst*136, 3707–3712 (2011). [PubMed: 21796289]
32. Park S, Zhang Y, Lin S, Wang T-H, Yang S, Advances in microfluidic PCR for point-of-care infectious disease diagnostics. *Biotechnol. Adv*29, 830–839 (2011). [PubMed: 21741465]
33. Melendez JH, Hardick J, Barnes M, Barnes P, Geddes CD, Gaydos CA, Molecular characterization of markers associated with antimicrobial resistance in *Neisseria gonorrhoeae* identified from residual clinical samples. *Sex. Transm. Dis*45, 312–315 (2018). [PubMed: 29465687]
34. Centers for Disease Control and Prevention (CDC), Update to CDC’s Sexually Transmitted Diseases Treatment Guidelines, 2006: Fluoroquinolones no longer recommended for treatment of gonococcal infections. *MMWR Morb. Mortal. Wkly Rep*56, 332–336 (2007). [PubMed: 17431378]
35. Kirkcaldy RD, Harvey A, Papp JR, del Rio C, Soge OO, Holmes KK, Hook EW, Kubin G, Riedel S, Zenilman J, Pettus K, Sanders T, Sharpe S, Torrone E, *Neisseria gonorrhoeae* antimicrobial susceptibility surveillance—The gonococcal isolate surveillance project, 27 Sites, United States, 2014. *MMWR Surveill. Summ*65, 1–19 (2016).
36. Allan-Blitz L-T, Hemarajata P, Humphries RM, Kimble M, Elias S, Klausner JD, Ciprofloxacin may be efficacious in treating wild-type gyrase a genotype *Neisseria gonorrhoeae* infections. *Sex. Transm. Dis*45, e18 (2018). [PubMed: 29465685]
37. Allan-Blitz L-T, Humphries RM, Hemarajata P, Bhatti A, Pandori MW, Siedner MJ, Klausner JD, Implementation of a rapid genotypic assay to promote targeted ciprofloxacin therapy of *Neisseria gonorrhoeae* in a large health system. *Clin. Infect. Dis*64, 1268–1270 (2016).
38. Melendez JH, Hsieh Y-H, Barnes M, Hardick J, Gilliams EA, Gaydos CA, Can ciprofloxacin be used for precision treatment of gonorrhea in public STD clinics? Assessment of ciprofloxacin susceptibility and an opportunity for point-of-care testing. *Pathogens*8, 189 (2019).
39. Cyr SS, Barbee L, Workowski KA, Bachmann LH, Pham C, Schlanger K, Torrone E, Weinstock H, Kersh EN, Thorpe P, Update to CDC’s treatment guidelines for gonococcal infection, 2020. *MMWR Morb. Mortal. Wkly Rep*69, 1911–1916 (2020). [PubMed: 33332296]
40. Goire N, Nissen MD, LeCornec GM, Sloots TP, Whiley DM, A duplex *Neisseria gonorrhoeae* real-time polymerase chain reaction assay targeting the gonococcal *porA* pseudogene and multicopy *opa* genes. *Diagn. Microbiol. Infect. Dis*61, 6–12 (2008). [PubMed: 18248938]
41. Geraats-Peters CWM, Brouwers M, Schneeberger PM, van der Zanden AGM, Bruisten SM, Weers-Pothoff G, Boel CHE, van den Brule AJC, Harmsen HG, Hermans MHA, Specific and sensitive detection of *Neisseria gonorrhoeae* in clinical specimens by real-time PCR. *J. Clin. Microbiol*43, 5653–5659 (2005). [PubMed: 16272500]
42. Verma R, Sood S, Bala M, Mahjan N, Kapil A, Sharma VK, Pandey RM, Samantaray JC, Evaluation of an *opa* gene-based nucleic acid amplification test for detection of *Neisseria gonorrhoeae* in urogenital samples in North India. *Epidemiol. Infect*140, 2110–2116 (2012). [PubMed: 22244259]

43. Maze MJ, Young S, Creighton J, Anderson T, Werno A, Nucleic acid amplification of the *opa* gene for detection of *Neisseria gonorrhoeae*: Experience from a diagnostic laboratory. *J. Clin. Microbiol*49, 1128–1129 (2011). [PubMed: 21248091]
44. Buckley C, Trembizki E, Donovan B, Chen M, Freeman K, Guy R, Kundu R, Lahra MM, Regan DG, Smith H, Whiley DM; GRAND Study Investigators, A real-time PCR assay for direct characterization of the *Neisseria gonorrhoeae* GyrA 91 locus associated with ciprofloxacin susceptibility. *J. Antimicrob. Chemother*71, 353–356 (2016). [PubMed: 26538505]
45. Gaydos CA, Cartwright CP, Colaninno P, Welsch J, Holden J, Ho SY, Webb EM, Anderson C, Bertuzis R, Zhang L, Miller T, Leckie G, Abravaya K, Robinson J, Performance of the abbot realtime CT/NG for detection of *Chlamydia trachomatis* and *Neisseria gonorrhoeae*. *J. Clin. Microbiol*48, 3236–3243 (2010). [PubMed: 20668135]
46. Schoepp NG, Schlappi TS, Curtis MS, Butkovich SS, Miller S, Humphries RM, Ismagilov RF, Rapid pathogen-specific phenotypic antibiotic susceptibility testing using digital LAMP quantification in clinical samples. *Sci. Transl. Med*9, eaal3693 (2017). [PubMed: 28978750]
47. Snodgrass R, Gardner A, Semeere A, Koppaarthi VL, Duru J, Maurer T, Martin J, Cesarman E, Erickson D, A portable device for nucleic acid quantification powered by sunlight, a flame or electricity. *Nat. Biomed. Eng*2, 657–665 (2018). [PubMed: 30906647]
48. Gootenberg JS, Abudayyeh OO, Kellner MJ, Joung J, Collins JJ, Zhang F, Multiplexed and portable nucleic acid detection platform with Cas13, Cas12a and Csm6. *Science*360, 439–444 (2018). [PubMed: 29449508]
49. Melendez J, Hardick J, Barnes M, Page K, Gaydos C, Antimicrobial susceptibility of *Neisseria gonorrhoeae* isolates in Baltimore, Maryland, 2016: The importance of sentinel surveillance in the Era of multi-drug-resistant gonorrhoea. *Antibiotics*7, 77 (2018).




Fig. 1. PROMPT platform design and operation.

(A) Comparison of standard of care for *N. gonorrhoeae* diagnosis (top) versus proposed clinical workflow with the PROMPT platform (bottom). (B) Schematic shows processing of patient swab samples. Steps include swab elution, mixing eluate with magnetic beads, and injection into the assay cartridge, followed by pH-mediated automation of nucleic acid binding, purification, and elution using transfer of the beads through cartridge wells containing preloaded buffers. (C) Photographs of the PROMPT assay cartridge. Cartridges are sealed before use with an adhesive tape, which is removed before injection of the sample into the first well. The instrument's faceplate detaches for loading the cartridge onto the PCR heat block. Remounting the faceplate with magnetic clasps aligns the cartridge with the instrument's magnet arm and fluorescence detector. (D) Photograph and schematic showing the instrument interfacing with a smartphone or tablet using Bluetooth for real-time reporting of PCR amplification. PCR was used to determine whether *N. gonorrhoeae* was

present in the clinical sample by *opa* amplification in the Cy5 channel (red) and whether the bacterial strain was susceptible to ciprofloxacin by wild-type (WT) *gyrA* amplification in the 6-carboxyfluorescein channel (blue). A.U., arbitrary units.

Author Manuscript

Author Manuscript

Author Manuscript

Author Manuscript




Fig. 2. Rapid PCR assay and instrumentation.

(A) Schematic of the internal components overlaid on a photograph of the PROMPT instrument. The assay cartridge mounts directly onto an aluminum PCR heat block, which is aligned with a fluorescence detector and servo-actuated magnet arm in the instrument. (B) Photograph of the magnet arm and schematic with cross-sectional view of the cartridge showing acrylic [polymethyl methacrylate (PMMA)], polytetrafluoroethylene (PTFE), and polypropylene (PP) plastic layers. The arm rotates and extends/retracts two opposing permanent magnets to transfer the magnetic beads between cartridge reagents stored in thermoformed wells. Red arrows indicate movement of the magnet arm, and blue arrows show the corresponding direction of bead transfer. (C) Comparison of traditional PCR and PROMPT PCR kinetics, with elevated primer and enzyme concentrations enabling less annealing and extension time of amplified products using the PROMPT platform. (D) Optimization of enzyme and primer conditions by cycle threshold (C_t) for PCR ($n = 1$ for each condition) with 1-s extension hold times using the *opa* assay. (E) Melt-curve analysis of *opa* amplicons using the 2 μ M primer conditions. (F) Schematic shows the miniaturized aluminum heat block, which decreases thermocycling time compared to traditional benchtop thermocycler routines. (G) Melt-curve analysis of PCR amplicons from a traditional benchtop PCR routine is compared to that for the PROMPT cartridge rapid PCR assay ($n = 1$ for each condition).




Fig. 3. PROMPT cartridge assay evaluation.

(A) Overall sample-to-answer workflow includes <1 min of hands-on time with DNA purification and with PCR amplification automated for a total time of 15 min. (B) Rapid thermocycling in the PROMPT instrument completes 40 cycles of PCR within 12 min during which a traditional thermocycling routine completes <10 cycles. (C) Cost breakdown for cartridge components. (D and E) Cartridge assay standard curves are shown using 10-fold serial dilutions run in triplicate ($n = 3$) from 1 to 10^5 CFU of *N. gonorrhoeae* input, with duplexed detection of *opa* (red) and wild-type *gyrA* (blue). (F) Time to result (minutes) using a live-detection algorithm versus *N. gonorrhoeae* CFU input to the cartridge. Error bars represent SD. (G) Cycle threshold (C_t) values for assay validation using reference strains of *N. gonorrhoeae* (ATCC, WHO), BCHD clinical isolates, related *Neisseria* species, and other organisms known to infect the urogenital tract. *N. gonorrhoeae* strains are color-coded by susceptibility to ciprofloxacin. No amplification is indicated by a C_t equal to 0.




Fig. 4. Clinical validation of the PROMPT platform.

(**A**) Results of *N. gonorrhoeae* detection by the PROMPT cartridge assay show *opa* amplification cycle threshold (C_t) for BCHD swabs ($n = 66$) and swabs collected in Kampala, Uganda ($n = 151$). Samples with no amplification in this plot are given a C_t set at 40 (gray, black). Horizontal dashed lines indicate cutoff C_t values used to determine positive amplification by the PROMPT cartridge assay. (**B**) Ciprofloxacin susceptibility of *N. gonorrhoeae* in BCHD samples was predicted by amplification of wild-type *gyrA* and is compared to culture results for antimicrobial susceptibility testing. (**C**) ROC curve for *N. gonorrhoeae* (NG) identification by *opa* amplification in all swab samples. (**D**) ROC curve for *N. gonorrhoeae* ciprofloxacin susceptibility by *gyrA* amplification. (**E**) Turnaround times for all swab samples using the complete 40-cycle PCR (blue) or determined by the live-reporting algorithm for *opa* gene amplification (red).